


Hardwood Leg Couch




MODEL 3620

- Dual frame construction with separate base & headrest sections
- Non-adjustable wedge headrest
- Solid, natural finish, hardwood legs
- Paper dispenser included
- 2" firm foam padding
- Heavyweight, knitbacked upholstery
- 250 lbs. load capacity under normal use

SPECIFICATIONS

3620-24	72" L x 24" W x 18" H
3620-27	72" L x 27" W x 18" H


California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Options



040 Paper Holder



016 Adjustable Wedge Headrest

Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black