

Small Up/Down Staircase

MODEL: 4-5500

- Small straight staircase for limited space
- Patient can ascend and descend without turning
- Available in 2 widths
- Requires a minimum 9' ceiling
- 500 lbs. bariatric load capacity under normal use



California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including vinyl chloride, which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

3

MODEL 4-5500-30 SPECIFICATIONS

Platform	30" W x 24" D x 12" H
Two Steps	30" W x 10" D x 6" H
Three Steps	30" W x 10" D x 4" H
Overall Size	34" W x 54" W x 47.75" H

MODEL 4-5500-36 SPECIFICATIONS

Platform	36" W x 24" W x 24" H
Two Steps	36" W x 10" D x 6" H
Three Steps	36" W x 10" D x 4" H
Overall Size	40" W x 54" D x 47.75" H