



# Commercial / Pharmaceutical Refrigerator

Conserv 24in COMMERCIAL/PHARMACEUTICAL Refrigerator 12.7cf WIFI 110V



## FEATURES

- No Frost
- USB Port
- Door Lock
- Adjustable Feet
- WiFi Compatible
- Interior LED Light
- Door Open Alarm
- Electronic Display
- Noise Level 43dB
- High-Low Temperature Alarm
- High-Low Temperature continuously displayed

## SPECIFICATIONS

Model	CMV 1300 G
Color	Blue on top/White
Capacity	12.7 cu.ft
Electrical	110V / 60Hz / 1.5A
Control Type	Buttons
Installation Type	Freestanding
Number of Shelves	7 (Grill)
Refrigerant	R600a
Temperature Control	Electronic
Temperature Range	35.6° F - 46.4° F
Weight (Net/Gross) lbs	157 / 194
Interior Dimensions HxWxD (inch)	57.3 x 19.7 x 23
Product Dimensions HxWxD (inch)	64 x 23.4 x 28
Package Dimensions HxWxD (inch)	65 x 25.2 x 30.1
Type of Shipping	LTL
Certification	ETL
Warranty	1 Year Parts & Labor

Keywords: Commercial Refrigerator, Pharmaceutical Refrigerator, Vaccine Refrigerator, Medical Refrigerator, Frost Free, Built-In Refrigerator



**Conserv Appliances** • 5006 Steffani Lane, Houston • TX 77041 • Toll Free 800-776-3538  
[www.ConservAppliances.com](http://www.ConservAppliances.com) • Email: [info@ConservAppliances.com](mailto:info@ConservAppliances.com)