




Commercial / Pharmaceutical Refrigerator



FEATURES

- No Frost
- USB Port
- Door Lock
- Adjustable Feet
- WiFi Compatible
- Interior LED Light
- Door Open Alarm
- Electronic Display
- High-Low Temperature Alarm
- High-Low Temperature continuously displayed

SPECIFICATIONS

Model	CMV 400
Color	Blue on top/White
Capacity	3.9 cu.ft
Electrical	110V / 60Hz
Control Type	Buttons
Installation Type	Freestanding/Built-In
Number of Shelves	3 (Grill)
Refrigerant	R600a
Temperature Control	Electronic
Temperature Range	35.6°F - 46.4°F
Weight (Net/Gross) lbs	82 / 89
Product Dimensions HxWxD (inch)	35 x 24 x 21.2
Package Dimensions HxWxD (inch)	35.8 x 25.2 x 23.2
Type of Shipping	FedEx / UPS
Certification	ETL
Warranty	1 Year Parts & Labor

Keywords: Commercial Refrigerator, Pharmaceutical Refrigerator, Vaccine Refrigerator, Medical Refrigerator, Frost Free, Built-In Refrigerator



Conserv Appliances • 10222 Georgibelle Drive, Suite 200, Houston, TX 77043-5249
Tel: 800-776-3538 • www.conservappliances.com • Email: info@ConservAppliances.com