



 **Convaid**

Convaid Cruiser

Lightweight, Compact-folding, 30° Fixed Tilt Wheelchair

Convaid Cruiser

The Convaid Cruiser is a lightweight, compact-folding, and highly durable wheelchair featuring Self-Tensio®, a built-in self-tensioning seating system that provides stable pelvic positioning. The Cruiser has a 30° fixed tilt, which improves postural control through upper body positioning for individuals with limited upper body strength and stability. The Cruiser is available in five sizes and provides years of growth through its seat depth extension tubes.




Lightweight, Compact-folding, and Highly Durable

The Convaid Cruiser is a 30° fixed tilt wheelchair that improves postural control through upper body positioning.



Tension Adjustable Strap Back

Allows for adjustment of the back according to the anatomical shape of the individual to provide optimal support.



Lightweight and Compact-folding

The folding mechanism allows the chair to fold into a tight and lightweight package.



Self-Tensio®

Provides key support at the pelvis. It provides a stable platform and distributes pressure throughout the contact area. It improves pelvic support and positioning without added weight and complexity.



5" of Seat Depth Growth

Handles the growth of the user for many years through extension tubes that can be added or changed easily without tools.

Standard Features

- 30° fixed tilt
- One-piece height adjustable push handle (CX10/12/14/16)
- Two-piece push handles (CX18 required)
- Tension adjustable back
- Cordura or Textilene upholstery
- Tool-less upholstery removal
- Tool-less adjustable seat depth growth of 5"
- Zippered pockets in seat & back upholstery (Cordura only: CX10/12/14/16)
- 2-point positioning belt
- Self-Tensio®
- Removable, height adjustable swing-away footplates
- 7.5" X 2" front, 11.5" X 2.5" Rear solid tires (CX10/12/14/16)
- 8" front and 11" rear solid tires (CX18 required)
- Quick release Rolko rear wheels (CX10/12/14/16)
- Foot operated wheel lock (CX10/12/14/16)
- Toggle wheel locks (CX18 required)
- Silver powder coated frame
- Closure strap
- Made in the USA

Options

Sizes

(Seat width)

- 10"
- 12"
- 14"
- 16"
- 18"

Upholstery Types


- Textilene
- Cordura
- Heavy duty Textilene
- Heavy duty Cordura

Tires


- 7.5" X 2" front, 12.5" rear solid or pneumatic knobby tires
- 12.5" rear solid knobby tires
- 8" X 2" front ribbed, 12.5" rear solid or pneumatic knobby tires

Upholstery Colors


Textilene




Wild raspberry



Teal green



Navy blue



Violet purple

Cordura



Princess pink



Cherry red



Apple green



Forest green



Aqua blue




Royal blue



Sassy purple



Chocolate brown



Panther black

Frame Options



Cruiser Size 18

*Manufacturer installation required



Cruiser Scout



Cruiser Planar

Accessories



**Extended headrest cover -
with and without windows
(canopy)**




Headrest cover (canopy)



Rain cover
Canopy required



Mosquito net
Canopy required



**Headrest extension for
Textilene upholstery**



**Headrest extension for
Cordura upholstery**



Occi-headwings



Padded headwings



Seat reducer insert



Full torso support vest



**Soft adjustable lateral
support single flap with
scoli strap**



**Soft adjustable lateral
support double flap with
scoli strap**

Accessories



H-harness with padded covers
3-point positioning belt required



Support cushion (general use)*¹



Position cushion (medial thigh support/anti thrust)*¹



Align cushion (lateral pelvis/lateral thigh support)*¹



3-point positioning belt




Lateral thigh support (adductor)



Medial thigh support (abductor)




Flip-up armrests
Available on the Cruiser Scout frame only.



Anatomic back*
Heavy duty reinforced upholstery required



Heavy-duty reinforced upholstery - Cordura or Textilene*



Attendant hand brake*



Upper extremity support surface (clear tray)

*Manufacturer installation required

¹Incontinence cushion liner available

Accessories



Under seat storage basket



Saddle bags




Wheelchair travel bag



Calf panel



Angle adjustable
footplates



Footplate securement
strap



Foot positioners



Castor Locks
(for 2" wide front tires)*



5-point harness with
padded covers
(CX10/12 only)*
Crash tested, WC-19
compliant



Q'straint transit lap belt



Sure-lok transit lap belt



Transit anchors*

*Manufacturer installation required

Specifications

Cruiser measurements		CX 10	CX 12	CX 14	CX 16	CX 18
	Seat width*	10"	12"	14"	16"	18"
	Seat depth	6-11"	8-13"	10-15"	11-16"	14-21"
	Seat back height^	21"(C) / 20"(T)	22"(C) / 20"(T)	25"(C) / 23"(T)	27"(C) / 24"(T)	25"(C) / 25"(T)
	Seat to footplate*	5-21"	5-25"	4-22"	4-22"	8-23"
	Seat to floor	23"	25"	21"	23"	26"
	Seat to back angle	85°/90°/95°	85°/90°/95°	85°/90°/95°	85°/90°/95°	85°/90°/95°
	Seat angle	30°	30°	30°	30°	30°
	Overall height	37"	38"	37"	40"	38"
	Overall length	48"	44"	49"	52"	55"
	Overall width	22"	23.5"	25.5"	27.5"	28.5"
	Folded height [◊]	17"	17"	17"	17"	17"
	Folded length	47"	48"	47"	53"	50"
	Folded width	15"	15"	16"	16"	15"
	Headrest extension height	6"	8"	8"	8"	10"
	Shoulder strap height [◊]	11-19"	12-20"	15-23"	17-25"	16-22"
	Push handle adjustments^^	28 - 43"	29 - 46"	28 - 45"	31 - 48"	N/A
	Weight capacity regular/transit	75/66 lbs	75/66 lbs	100/100 lbs	170/170 lbs	250/170 lbs
	Weight of chair without front rigging	27 lbs	27 lbs	28 lbs	32 lbs	30 lbs

* Plus or minus .5 inches due to fabric
 ◊ Models with Textilene fabric, shoulder strap height is 2 inches shorter
 ^ Seat back height
 C - Cordura upholstery
 T - Textilene upholstery
 ^^ Adjustable push handle is not available on 18 inch models

BR1900-ENG-Rev.00 2022-07-19 | Reserving our rights for printer's error discontinued products. Please find an updated edition on our website © 2022 Etac All rights reserved

For the latest news and continuously updated
product information visit: www.etac.us.com

Etac
+1 844-US-MOBILITY (844-876-6245)
convaidsales.us@etac.com
www.etac.us.com

DCN0956.0


Creating Possibilities