

PreserveTech Aquachair Transfer Accessory with Curtain Guard Protection

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

Enhance your Aquachair Bathing System with the bench attachment that makes getting out of the tub easier and safer. The PreserveTech Aquachair Transfer Accessory with Curtain Guard allows you to easily slide over the side of the tub without standing. For added protection against slips and falls, its "U" shape cut-out design keeps the shower curtain liner inside the tub to prevent water from spilling during bathing.

- Enhance Your Aquachair with Safer Entry &
- Easily installs to the PreserveTech Aquachair Bathing System with Bidet in minutes
- Bench design allows user to slide over side of tub from a safer seated position
- Patented "U" shape cut-out design keeps shower curtain liner in the tub and prevents water spillage
- Adjust leg height with the push of a button
- Nonslip feet offer added security on bathroom floors

PRODUCTS	Description
RTL12A004-TR	AQ Transfer Bench Accessory, 1EA, 1CV
SPECIFICATIONS	Description
Warranty	1yr Limited
Product Weight	3.85 lb
Material	PP+Aluminium
Weight Capacity	350 lb
Color	White
Dimensions Assembled	13.43" (W) X 12.28" (L) X 18.58" (H)
Retail Box Dimensions	17.32" (L) X 14.37" (H)
Retail Box Weight	2.02 lb
CA Prop 65	no
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]