

Rear Glide Brakes with 3" Wheels

HCPCS: A9900 OR E0155



FIGURE A



- For use on back of wheeled walker
- 3" wheel glides along surface until weight is applied to hand grip, which activates rubber-tip brake
- Allows for 8 height adjustments
- Weight Capacity 350 lb

PRODUCTS

Description

10106	Glide Brakes w/ 3in Wheel, 1 pr/bx, RTL
-------	---

SPECIFICATIONS

Description

Weight Capacity	350 lb
-----------------	--------