

Swivel Wheel with Lock & Two Sets of Rear Glides

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Can be locked for fixed wheel or unlocked for a swivel wheel (Figure A - B)
- Allows for 5 adjustments
- Rubber wheels allow walker to roll easily and smoothly over irregular surfaces
- Comes with rear glide caps, (item #10107) and glide covers (item #10107C), allowing use on all surfaces
- Weight Capacity 350 lb

PRODUCTS

Description

10115	Swivel-Lock Wheel,5",1 pr/bx, RTL
-------	-----------------------------------

SPECIFICATIONS

Description

Weight Capacity	350 lb
-----------------	--------