

5" Bariatric Walker Wheels with Two Sets of Rear Glides

HCPCS: E0149



FIGURE A



FIGURE B



FIGURE C

- Converts heavy-duty walker into wheeled walker
- Allows for 8 height adjustments
- 10118T Adds 4"
- Dual rubber wheels allow walker to roll easily and smoothly over irregular surfaces
- Comes with rear glide caps (item #10107) (Figure A) and glide covers (item #10107C) (Figure B), allowing use on all surfaces
- For use with items # 10220-2 and 10223-2

PRODUCTS	Description
10118SV	Walker Wheels,HD,Silver,1PR/bx,RTL
SPECIFICATIONS	Description
Weight Capacity	500 lb