

Bariatric Aluminum Folding Walker- Two Button

HCPCS:



- For Bariatric use
- Steel legs and side braces
- Precision design provides additional strength while adding minimal weight
- Wider and deeper frame design accommodates individuals up to 600 lb
- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Vinyl-contoured hand grip

PRODUCTS	Description
10220-1	Walker,Oversized,Adt.,1/ea

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Handle Height	32" - 39"
Weight Capacity	600 lb
Carton Shipping Weight	12 lb
Overall Width	23"
Weight	9 lb