

Bariatric Steel Forearm Crutch

HCPCS:



- Vinyl-coated, ergonomically contoured arm cuffs molded for comfort and stability
- Leg and forearm sections adjust independently for optimal sizing
- Vinyl hand grips are comfortable and durable
- Extra-large tips for added stability

PRODUCTS	Description
10403HD	ForearmCrutch,Bariatric,Adult,1/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	500 lbs
Weight	5.7 lbs.
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]