

Revolving, Adjustable-Height Stool

HCPCS:



- Smooth screw height adjustment from 17.5"24"
- Solid steel screw shaft
- 1" chrome plated steel tubular construction
- 14" seat diameter
- 4" thick padded seat (black)
- Four hooded casters
- Optional seat colors available
- Minimum quantity required

PRODUCTS

Description

13034	Stool, Round, Blk,Wheeled,1/cs
-------	--------------------------------

SPECIFICATIONS

Description

Warranty	Limited Lifetime
----------	------------------

Weight Capacity	250 lb
-----------------	--------

Height Adjustment	17.5" - 24"
-------------------	-------------