

Mayo-Instrument Stand

HCPCS:



- Single post stand with "U" base design provides support and versatility
- Removable stainless steel tray measures 19" x 13"
- Tray height adjusts from 32" - 49" with lock
- Two, 3" casters provide additional mobility

PRODUCTS

Description

13035	Mayo-Inst.Stand,Single Post,1/cs
-------	----------------------------------

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Height	32" - 49"
Base Dimensions	22" (W) x 18.5" (D)
Carton Shipping Weight	11 lb
Weight	9 lb

ACCESSORIES

Description

13035T	Tray for 13035, 13045,13071,1/ea, EA
--------	--------------------------------------