

## Mayo-Instrument Stand

HCPCS:



- Double post provides added support for use in surgical procedures
- Removable stainless steel tray measures 19" x 12.5"
- Adjustable height from 32.5" - 52"
- Two hooded casters

| PRODUCTS               | Description                           |
|------------------------|---------------------------------------|
| 13045                  | Mayo-Instrument Stand,DoublePost,1/cs |
| SPECIFICATIONS         | Description                           |
| Warranty               | Limited Lifetime                      |
| Height                 | 32.5" - 52"                           |
| Base Dimensions        | 21" (W) x 19" (D)                     |
| Carton Shipping Weight | 14 lb                                 |
| Weight                 | 11 lb                                 |
| ACCESSORIES            | Description                           |
| 13035T                 | Tray for 13035, 13045,13071,1/ea, EA  |