

Free-Standing Trapeze

HCPCS: E0912



- Designed to assist larger individuals in changing positions while in bed and aid in transferring in and out of bed
- Easy-to-assemble tool-free design
- Lightweight design for easy delivery
- Chain is adjustable from edge of headboard to the end of boom by sliding to desired position

PRODUCTS	Description
13049SV	Bariatric Trapeze Silver Vein,650lb,1/ea
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	650 lb
Boom Length	37"
Carton Shipping Weight	83 lb
Weight	74 lb