

Mayo-Instrument Stand

HCPCS:



- Casters provide additional mobility and convenience
- Removable stainless steel tray measures 19" x 13"
- Tray height adjusts from 29.5" - 47" with lock

PRODUCTS	Description
13071	Instrument Stand,Mobile,5caster,1/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Height	29.5" - 47"
Base Dimensions	21" (W) x 21" (D)
Carton Shipping Weight	12 lb
Weight	10 lb
ACCESSORIES	Description
13035T	Tray for 13035, 13045,13071,1/ea, EA