

Battery-Powered Patient Lift

HCPCS: E0635



- Heavy-gauge construction
- Jumbo actuator pump with emergency button can switch to manual mode to lower patient down safely
- Fast, easy charge from AC outlet Lift will not operate while plugged in For home-care use
- Battery does not need to be removed to be charged
- Audible warning when battery power falls below 50%

PRODUCTS	Description
13240	Patient Lift - Electric, 1/cs

SPECIFICATIONS	Description
Warranty	Limited Lifetime on frame 1 year on components and electrical parts
Weight Capacity	450 lb
1 Year Electronics	Battery, Hand control, Actuator
Base Clearance	2"
Boom Height	27.5" - 76"
Carton Shipping Weight	122 lb
Casters	5" Rear
Motor	Linak Motor
Weight	108 lb