

Bariatric Battery-Powered Lift

HCPCS: E0635



- To charge, simply plug in unit Battery does not need to be removed to be charged
- Lift will not operate while plugged in Ideal for home-care use
- Jumbo actuator pump with emergency button can switch to manual mode to lower patient down safely
- Base legs adjust easily and lock securely into open position with padded shifter handle
- Extra wide coated cradle with 360° rotation with 4-point hook up (Can fit under bed)
- Weight Capacity 600 lbs

PRODUCTS	Description
13244	Bariatric Lift,Elect,Homecare Style,1/cs

SPECIFICATIONS	Description
Warranty	Limited Lifetime on frame 1 year on components and electrical parts
Weight Capacity	600 lbs
Base Clearance	2"
Boom Height	32" - 75.5"
Carton Shipping Weight	148 lbs.
Casters	3" Front 4" Rear
Motor	TiMotion Motor
Weight	132 lbs.