

Pull Cord Alarm

HCPCS:



- Economical way to notify staff that a patient has moved
- Alarm easily secures to a bed or chair with a self-contained clip
- Activation cord is adjustable from 18"36"
- When user gets up, the pin connector pulls from the alarm unit activating the signal
- Deactivate with On/Off switch
- Two volume settings: 97 dB and 103 dB
- Requires one 9V battery (Included)

PRODUCTS	Description
13602	Cord Alarm,w/Pin,1/ea
SPECIFICATIONS	Description
Warranty	Limited 1 Year
ACCESSORIES	Description
136DC	ALARM AC ADAPTER, 1/EA, 25/CV, CS