

Pressure-Sensitive Chair and Bed Patient Alarm

HCPCS:



FIGURE A

- Alerts care giver with audio alarm when patient gets out of a chair or bed
- Pressure-sensitive pad connects to the audio alarm
- The alarm comes with a safety alert, On/Off switch and low battery warning
- Comes with 9V Battery (Included)

PRODUCTS	Description
13605	Patient Alarm,Chair,1/ea
13606	Patient Alarm,Bed,1/ea
13608	Chair Size Patient Alarm w/Res,1/ea
13609	Bed Size Patient Alarm w/Reset,1/ea
SPECIFICATIONS	Description
Warranty	1 Year (alarm) 6 Months (sensor pad)
ACCESSORIES	Description
136DC	ALARM AC ADAPTER, 1/EA, 25/CV, CS