

Standard Telescoping Full-Length Side Rail, Brown-Vein Finish

HCPCS: E0310





FIGURE A



FIGURE B




FIGURE C



FIGURE D

- Attractive brown-vein finish (Figure A - B)
- Constructed of 1" steel
- Telescopes with movement of bed
- Easily adjusts up or down with spring-loaded release
- Spring-loaded crossbars are easy-to-install and remove (Figure C)

PRODUCTS	Description
15001ABV	Bed Rails,Full,BrownVein,4 Bar,1/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]
ACCESSORIES	Description
15050BV	Cross Bars, Universal,Brown Vein, 1/pr, CS