

Bariatric Foam Mattresses

HCPCS: E0184-1,E0184-1




FIGURE A

- Nylon top cover is fluid-resistant and easy to clean for extended life and service
- Vinyl bottom is waterproof
- Fire Retardant: Mattress meets the requirements of 16 CFR Part 1632 Standard for the Flammability of Mattresses and Mattress Pads and 16 CFR Part 1633 (Federal Flammability (open flame) standard for mattress sets) when used without a foundation
- High-Density foam mattress for the bariatric patient
- Solid foam inner core provides comfort, support and durability in a one-piece construction
- Latex free

PRODUCTS	Description
15301	Bariatric Mattress, 42in,1/ea
15301-84	Bariatric Foam Matt,84"x42"x6.5",1/ea
15310-84	Bariatric Foam Matt,84"x48"x6.5",1/ea
SPECIFICATIONS	Description
Warranty	1 Year
Product Dimensions	15301: 80(L) x 42(W) x 6.5(H) 15301-84: 84(L) x 42(W) x 6.5(H) 15310-84: 84(L) x 48(W) x 6.5(H)
Weight Capacity	15301: 585 lbs 15301-84: 585 lbs 15310-84: 665 lbs
FireRating	CFR 16 Part 1633 Federal Fire Code
Side Rails	No
Cover Material	Fluid-Resistant Nylon Top Waterpoof Vinyl Bottom
Limited Warranty	1 Year
Mattress Layers of Foam	1
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]
ACCESSORIES	Description
15301-84CVR	Cover Bari Foam Matt,84"x42"x6.5",1/ea, EA
15310-84CVR	Cover Bari Foam Matt,84"x48"x6.5",1/ea, EA
15312-84CVR	Cover Bari Foam Matt,84"x54"x6.5",1/ea, EA
15301C	Bariatric Mattress Cover, 42", 1/ea, EA
15310C	Mattress Cover 48" 1/cv, CS