

Heavy Duty 50 PSI Compressor

HCPCS: E0565



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Newly designed compressor vibration mounts, clamps and lock washers, and wire harnesses improve upon an already durable design
- Ultra-powerful compressor for both humidification and high flow, high-pressure nebulization of viscous medications
- Designed to provide efficient, continuous high-pressure performance
- Lightweight and easy-to-clean

PRODUCTS

Description

18450	Heavy Duty 50PSI Nebulizer,1/ea
-------	---------------------------------

SPECIFICATIONS

Description

Warranty	Limited 1 Year
Flow Rate	11~41 LPM
Power Consumption	160 W

ACCESSORIES

Description

MASK 007A	Adult Tracheotomy Mask 50/case, CS
MASK 007P	PediTracheotomy Mask 50/CS, CS
MASK 001P	Pediatric Aerosol mask CS/50, CS
MASK 001A	Adult Aerosol mask CS/50, CS