

MD by Drive Air-Go Walker

HCPCS:




FIGURE A




FIGURE B



FIGURE C

Differentiate your offerings with the Air-Go Walker--delivering comfort, ease, and freedom beyond traditional walkers. With oversized, cushioned Airgonomic handles, water-dispersing treaded wheels, and a padded mesh seat, it gives users both freedom of movement and confidence in every step.

- **Airgonomic Comfort:** Oversized hand grips are easy to grasp and feature integrated air bladders that cushion the hands.
- **Convenient Seating:** Integrated, open-and-go mesh seat allows for rest anytime.
- **Natural Stride:** Curved $\frac{3}{4}$ -length seat design allows users to walk comfortably without interference.
- **Easy, Fast Fold:** The simple scissor fold design allows for compact storage and easy transportation.
- **Improved Traction & Stability:** Large, water-dispersing treaded front wheels grip the earth help maneuver through diverse terrain.
- **Fully Assembled Out of the Box:** Ready to use immediately, requiring no tools or setup for added convenience.

PRODUCTS	Description
MD30100BK	Air-Go Walker, Onyx
MD30100BL	Air-Go Walker, Sapphire
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	8 lb
Material	Anodized Aluminum
Seat Dimensions	21.5" (W) x 8.8" (D)
Handle Height	33.1" to 36"
Weight Capacity	300 lb
Height	Offering adjustments for users 4'9" - 6'7"
Seat To Floor Height	20 - 23
Dimensions Assembled	25.6 (W) x 25.2 - 26.9 (L) x 33.1 - 36 (H)
Folded Dimensions	7.5" (W) x 26.2" (L) x 20.5" (H)
Retail Box Dimensions	21.5 (L) x 7.1 (W) x 28.5 (H)
Wheel Size	6 (D) x 1 (W)