

## Aluminum Crutches with Accessories

HCPCS: E0114,E0114



FIGURE A




FIGURE B




FIGURE C

- Lighter than wood, the standard aluminum crutch provides stability and durability
- Double extruded center tube provides additional strength to weight-bearing area
- Pushpin adjustment makes it easy-to-adjust crutch height in 1" increments
- Non-skid, jumbo-size vinyl contoured tips provide excellent traction
- Easy wing nut adjustments for hand grips
- Comfortable, durable accessories are latex free
- Weight Capacity 350 lb on aluminum crutch (excluding pediatric)
- Piggy back label supplied on product hangtag (excluding pediatric 10416-1)

| PRODUCTS        | Description                                                                                                                                                                                                                                                 |
|-----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 10400-8         | Crutch, Adult, 8pr/cs                                                                                                                                                                                                                                       |
| 10401-1         | Crutch,Alum.Youth,1pr/bx                                                                                                                                                                                                                                    |
| 10401-8         | Crutch, Youth, 8pr/cs                                                                                                                                                                                                                                       |
| 10402-8         | Crutch, Tall, 8pr/cs                                                                                                                                                                                                                                        |
| SPECIFICATIONS  | Description                                                                                                                                                                                                                                                 |
| Warranty        | Limited Lifetime                                                                                                                                                                                                                                            |
| Weight Capacity | 350 lb                                                                                                                                                                                                                                                      |
| CA Prop 65      | [WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <a href="http://www.P65Warnings.ca.gov">www.P65Warnings.ca.gov</a> .] |