

## Bariatric Platform Walker/Crutch Attachment

HCPCS: E0154



FIGURE A



- Designed for those who are unable to grip the walker
- Cradles forearm with soft vinyl padding for comfortable weight bearing
- Large top bracket allows the platform attachment to be placed over the hand grip of the walker
- Adjustable strap securely holds arm in place
- Attachment for Underarm Crutch (Figure A) for item #10105 must be ordered separately (Cannot be used with Bariatric model)
- Bariatric Weight Capacity 500 lb

| PRODUCTS  | Description                             |
|-----------|-----------------------------------------|
| 10105HD-2 | Walker Platform Attatcmt,HeavyDuty,2/cs |

| SPECIFICATIONS         | Description      |
|------------------------|------------------|
| Warranty               | Limited Lifetime |
| Weight Capacity        | 500 lb           |
| Height                 | 40.5" - 55"      |
| Carton Shipping Weight | 9 lb             |
| Weight                 | 4 lb             |