

Blue Streak Wheelchair

HCPCS: K0001-1,K0001-1



- Blue Powder-coated frame, black cross brace
- Flip-back, padded desk-length arm
- Nylon upholstery is durable, lightweight, attractive and easy to clean
- Solid Polyurethane tires mounted on composite wheels provide durability and low maintenance
- Comes with push-to-lock wheel locks
- Comes with calf strap (-SF models only)
- Comes with swing-away footrests or elevating leg rests

PRODUCTS	Description
BLS16FBD-ELR	BlueStreak16"D/L FlipBackArms,Elev,1/ea
BLS16FBD-SF	BlueStreak16"D/L FlipBackArms,SwgFt,1/ea
BLS18FBD-ELR	BlueStreak18"D/L FlipBackArms,Elev,1/cs
BLS18FBD-SF	BlueStreak18"D/L FlipBackArms,SwgFt,1/cs
BLS20FBD-ELR	BlueStreak 20"D/FlipBackArms,EivLeg,1/ea
BLS20FBD-SF	BlueStreak 20"D/FlipBackArms,SwgFt,1/ea

SPECIFICATIONS	Description
Warranty	Limited 3 Year
Product Weight	41.25 lb
Weight Capacity	300 lb
Seat To Floor Height	19.5"
Carton Shipping Weight	47 lb
Back Height	16"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]

ACCESSORIES	Description
STDS833N	Rear Anti Tipper w.wheels SSP&CS,1/pr, PR
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS801GT	WheelLockExtension,ViperGT,SSP2,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
WASL	Limb Support,Left,1/cs, CS

ACCESSORIES**Description**

WASR	Limb Support, Right,1/cs, CS
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS831	Heel Loops, 1pr, BX