

CarbonTouch™ Power Wheelchair

Travel-Ready Carbon Fiber Mobility



Meet the growing demand for travel-friendly mobility with the CarbonTouch Power Wheelchair.

With an ultra-lightweight carbon fiber frame and high weight capacity, this airline-ready power wheelchair brings premium comfort to more customers while maintaining exceptional portability for active lifestyles.

Product Specs

Drive Range	10.5 miles
Top Speed	3.70 mph
Turning Radius	27.6"
Ground Clearance (Middle)	2.4"
Width Between Arms	20"
Armrest Type	Flip-Back, Depth-Adjustable
Leg rest Type	Footplate
Seat Dimensions	17" (W) x 18" (D)
Seat to Floor Height	19"
Folded Dimensions	33.3" (H) X 23" (W) X 15.5" (D)

Benefits & Features

Ultra-Portable & Travel-Friendly

At just 36.4 pounds with a one-touch folding mechanism, and FAA-compliant lithium battery, the CarbonTouch is easy to carry and perfect for trains, planes, cruise ships, and more.

Accommodates More Users

A 330-pound weight capacity and 17" (W) x 18" (D) seat accommodate more patients than competing models. Easy transition between powered and manual modes allows flexible use across care settings.

Safety & Convenience Included

Anti-tippers, safety belt, and removable joystick with integrated USB charging port all come standard. Dual charging options allow for on-board and off board charging.

Durable Design

Sleek carbon fiber frame offers superior durability while being significantly lighter than traditional materials.



Battery Type (included)	Lithium Ion 24V 12Ah
Battery Charger Type	Off Board & Controller Port
Wheel Size	7" Front 8" Back
Base Weight	40.8 lb. 36.4 lb. without Battery
Battery Weight	4.4 lb.
Controls	Removable Joystick
Weight Capacity	330 lb.
Warranty	1 year motor, controller and CPU 5 years frame, 1 year battery, 6 months wearable parts

Item #	Description	UPC	UOM
CFPC17FS	CarbonTouch™ Power Wheelchair	822383030319	1/ea