

Seat2Go Positioning Seat

HCPCS: E1399-1,E1399-1



FIGURE A



- Contoured seat prevents child from sliding forward and improves seating posture
- 95° seat to back angle
- Double-strap trunk support and a positioning belt
- Two sets of straps secure Seat2Go to any chair
- When not in use, the seat straps can be stored under the seat behind the clear flap and the back straps in the back pocket
- Can be used as an insert for a single or double stroller
- Weighs 5 lb - 7 lb

PRODUCTS	Description
CRS 2000	Seat2Go, small,1/cs
CRS 3000	Seat2Go, Medium,1/cs
SPECIFICATIONS	Description
Warranty	Limited 1 Year
Product Weight	5 lb 7 lb
Carton Shipping Weight	8 lb 11 lb
Overall Width	8" 12"
Back Height	12" 14"
ACCESSORIES	Description
CRS 8025	Abductor, Seat2GO sm/lg, 1/bx, BX
CRS 8000	Headrest, Seat2Go & Back2Go, 1/ea, BX