

Cane Seat

HCPCS: E0100



FIGURE A



FIGURE B

- Provides a comfortable seat to rest on when open, and a sturdy support cane when closed (Figure A)
- Manufactured with sturdy, extruded aluminum tubing
- Tripod design with vinyl contoured tipped legs
- 9" diameter seat
- 33" tall when closed
- Weighs 2.5 lb
- Weight Capacity 250 lb

PRODUCTS	Description
RTL10365-ADJ	Adjustable Height Cane Seat,1 EA
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	250 lb
Seat To Floor Height	21"
Cane Height	34" - 38"