

Universal Cup Holder

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- New clamp allows cup holder to swivel to keep contents level
- Easily attaches to wheeled walkers, rollators, transport chairs/wheelchairs with tubing diameter from 7/8" to 1"
- Special clamp insert can be added or removed as necessary to accommodate different size tubing
- Made with durable, lightweight plastic

PRODUCTS

Description

STDS1040S	Cup Holder, Swivel, 1 BX, 6/CS
-----------	--------------------------------

SPECIFICATIONS

Description

Dimensions	3.5" (W) x 4" (L).
------------	--------------------