

Clinical Care Recliner

HCPCS:



- Ideal for use in dialysis, blood collection, respiratory and cardiac care
- Chair can be positioned in 4 positionsupright (Figure A), TV (Figure B), Full Recline (Figure C) and Trendelenburg
- Deep recline and Trendelenburg position can only be attained by an attendant
- Gas cylinder controls deep recline and Trendelenburg position
- Side panels "pop-off" for easy cleaning and maintenance
- Large, blow-molded side trays includes recess for cup
- 5" Casters (2 with locks) make it easy-to-move and secure recliner
- Gas cylinder release handle for quick back release
- Release handle can be operated from either side of the recliner
- Wide, comfortable armrests are angled to provide support in any reclined position
- Upholstered arm supports provide extra comfort in recline position
- Neutral-colored, sturdy steel frame is easy to clean and maintain
- Comfortable built-in headrest
- Foam and vinyl meet California Technical Bulletin 117 Standard

PRODUCTS	Description
D577-CHAR	Clinical Recliner, Charcoal
D577-CHOC	Clinical Recliner, Chocolate
D577-TAN	Clinical Recliner, Tan
SPECIFICATIONS	Description
Warranty	Limited Lifetime on Frame, 1 Year on Upholstery and Moving Parts
Product Weight	115 lbs.
Seat Dimensions	21.5" (H) X 21.5" (W) X 20" (D)
Weight Capacity	250 lbs
Carton Shipping Weight	135 lbs.
Width Between Arms	22"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov]