

Deluxe Fly-Weight Aluminum Transport Chair with Removable Casters

HCPCS: E1038



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Weighs only 19 lb
- Aluminum frame for lightweight strength and durability
- Easy folding and storage in closet or trunk
- Composite wheels are lightweight and maintenance free
- 8" wheels in front and rear easily remove without tools
- Deluxe back release
- Deluxe, all-aluminum rear wheel locks
- Comes with carry pocket on backrest
- Tool-free, swing-away footrests included
- Seat belt included for added safety
- Flared arms allow the back to be folded virtually flat for easy storage and transport
- Soft, plush nylon upholstery
- Convenient cup holder can easily attach to chair (sold separately Model# STDS1040S)
- Weight Capacity: 300 lb

PRODUCTS	Description
RTLFW19RW-RD	Flyweight,19in Rem Whls,Red,RTL,1 EA
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	19 lb
Product Dimensions	33" (L) w/ riggings x 22" (W) x 37.25" (H)
Seat Dimensions	19" (W) x 16" (D)
Weight Capacity	300 lbs
Seat To Floor Height	19"
Carton Shipping Weight	24 lbs
Folded Dimensions	33" (L) w/ riggings x 9" (W) x 37.25" (H)
Back Height	18"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]
ACCESSORIES	Description
STDS1040S	Cup Holder, Swivel,1BX,6/CS, BX