

## Deluxe Folding Walker, Two Button

HCPCS:



FIGURE A



- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Unique U-shaped frame design creates greater clearance
- The patient is able to bring the walker closer to assist in standing because of the improved crossframe (Figure A)
- Sturdy 1" diameter, aluminum construction ensures maximum strength while remaining lightweight

PRODUCTS	Description
10200-1	Walker,2Btn.Adult,1/ea
10200-4	Deluxe Two Button Walker, 4/CS
10201-4	Deluxe Two Button Junior, 4/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	350 lb
Carton Shipping Weight	26   24 lb
Overall Width	20"
Weight	5.5 lb
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <a href="http://www.P65Warnings.ca.gov">www.P65Warnings.ca.gov</a> .]