

Deluxe, Trigger Release Folding Walker

HCPCS:



FIGURE A



- Trigger release feature for opening and closing allows users to close the walker without removing their hands from the handle and is ideal for individuals with limited finger dexterity (Figure A) Trigger release can be reached with hands of all sizes
- Durable, Composite Trigger Release resists cracking
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight
- Patented

| PRODUCTS | Description |
|------------------------|----------------------------------|
| 10224-4 | Walker,Dual,TriggerRel,Adlt,4/cs |
| 10225-4 | Walker,DualTriggerRel.Jun.,4/cs |
| SPECIFICATIONS | Description |
| Warranty | Limited Lifetime |
| Handle Height | 32" - 39" 25" - 32" |
| Weight Capacity | 350 lb |
| Carton Shipping Weight | 25 lb |
| Overall Width | 23.5" |
| Weight | 6 lb |