

First Class School Chair

HCPCS: E1399-1,E1399-1



FIGURE A




FIGURE B

- Seat height and depth adjustable solid seat with a 15° anterior or posterior tilt
- Armrests provide lateral limitations as well as a stable resting surface for arms
- Flip-up armrests allows for easy transfers and are height adjustable
- The tray accessory mounts onto the armrest
- Padded Hip Belt 2-point hip belt provides pelvic support and maintains proper hip positioning

PRODUCTS	Description
FC 2000N	FirstClass SchoolChair,Sm,NewStyle,1/ea
FC 4000N	FirstClass SchoolChair,Lg,NewStyle,1/ea
SPECIFICATIONS	Description
Product Weight	13 lbs. 18 lbs.
Weight Capacity	100 lbs 200 lbs
Seat To Floor Height	10-14" 15 - 23"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]
ACCESSORIES	Description
FC 2028	Footrest, FC2000, 1 pr/bx, BX
FC 4028	Footrest, FC4000,1/pr/bx, BX
FC 2024	Dining Tray,FC2000,1/bx, BX
FC 4024	Dining Tray, FC 4000,1/bx, BX
FC 8000N	Multi-Axis Headrest, 1/ea, EA
FC 2030N	Anti Tippers, FC2000N, 1 pr, PR
FC 2029N	Legs w/Casters,FC2000,1 set/bx, BX
FC 4029N	Legs w/Casters,FC4000,1 set/bx, BX
FC 4030N	Anti Tippers, FC 4000N, 1 pr, PR
FC 8027N	Hip Guide, NewStyle, FC, 1 pr, PR
FC 8025N	Support Kit, Sm, FC2000N, 1 set/bx, BX
FC 8025 LN	Support Kit, Lg, FC4000N, 1 set/bx, BX