

Folding Quad Cane, Small Base, Charcoal, 1EA, 12/CV

HCPCS: E0105,E0105



FIGURE A




FIGURE B



FIGURE C



FIGURE D

- The Folding Quad Cane features a small footprint and folding design, making it ideal for display, storage and transport. Equipped with a four-point base and offset handle, the cane provides added stability and adjusts from 30-39 to accommodate varying user heights. The cane is made of extruded aluminum for added strength, while the contoured handle offers a comfortable grip.
- Folds and disassembles for small footprint and easy storage
- Made of sturdy, extruded aluminum tubing
- Four-point base and offset handle for additional stability and support
- Contoured handle for comfortable grip
- Handle height adjusts with the push of a button from 30 to 39
- Locking ring eliminates rattling and provides added safety

PRODUCTS

Description

RTL10310FCH	Folding Quad Cane, Charcoal, 1EA, 12/CV
-------------	---