

Folding Quad Cane, Small Base

The Folding Quad Cane features a small footprint and folding design, making it ideal for display, storage and transport. Equipped with a four-point base and offset handle, the cane provides added stability and adjusts from 30–39" to accommodate varying user heights. The cane is made of extruded aluminum for added strength, while the contoured handle offers a comfortable grip.



CONTOURED,
COMFORTABLE
HANDLE



FOLDS & DISASSEMBLES
FOR SMALL FOOTPRINT
AND EASY STORAGE



STURDY ALUMINUM
TUBING SUPPORTS
300 LB.



FOUR-POINT BASE
& OFFSET HANDLE
FOR STABILITY



HANDLE HEIGHT
ADJUSTS FROM
30–39"

RETAIL CARTON HAS
A 50% REDUCTION IN
FOOTPRINT FROM A
STANDARD QUAD CANE



ITEM #	PRODUCT DESCRIPTION	COLOR	UPC	UOM
RTL10310FCH	Folding Quad Cane, Small Base	Charcoal	822383008547	1ea/12cv

PRODUCT SPECS

Minimum/Maximum Height	30"/39"
Product Dimensions	5.75" (W) x 24" (D)
Product Weight	1.45 lb.
Color	Charcoal

BOX SPECS

Retail Box Dimensions	9" (W) x 3.5" (D) x 17" (H)
Retail Box Weight	1.45 lb.
Weight Capacity	300 lb.
Warranty	Limited Lifetime*

*Refer to instruction guide or website for warranty details.