

Folding Steel Commode

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- Easy, tool-free assembly
- Easily opens and folds for storage and transportation
- Durable plastic snap-on seat and lid
- Powder-coated steel-welded construction
- Plastic arms for added comfort
- Includes 7.5 qt commode bucket with carry handle, cover and splash shield
- Retail Packaged
- Made using recycled materials

PRODUCTS	Description
RTL11158KDR	Commode, Folding Frame, RTL,1 EA
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	13.75" (W) x 16.25" (D)
Weight Capacity	350 lb
Carton Shipping Weight	17.5 lb
Weight	13 lb
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]