

Steel Forearm Crutches

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Vinyl-coated, ergonomically contoured arm cuffs molded for comfort and stability
- Leg and forearm sections adjust independently for optimal sizing
- Vinyl hand grips are comfortable and durable
- Extra-large tips for added stability
- Ortho K Handle is hard plastic

PRODUCTS

Description

10403	Forearm Crutch, Adult, 1pr/bx
10405	Forearm Crutch, Tall 1pr/bx

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Weight Capacity	300 lb
Weight	2.5 - 2.6 lb
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]

ACCESSORIES

Description

10439-8	CrutchTips,10400,10406,10405Gry, 8/csRTL, CS
---------	--