

Moxie GT Gait Trainer

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Height adjustable in 1" increments
- Accessories mount tool free and can be adjusted according to the user's needs
- Can be used in anterior or posterior position
- Medium and Large frames fold for easier transportation
- Swivel casters can be locked in non-swivel position
- Directional lever allows for both forward and reverse mobility, or forward mobility only
- Wheel locks on all 4 casters to keep gait trainer stationary during transfers
- Small and Medium sizes available in blue and green
- Large size available in grey

PRODUCTS

Description

GT1000-2GB	Moxie GT Gait Trainer, Small, Blue
GT2000-2GB	Moxie GT Gait Trainer, Medium, Blue
GT2000-2GG	Moxie GT Gait Trainer, Medium, Green
GT3000-2GGY	Moxie GT Gait Trainer, Large, Grey

SPECIFICATIONS

Description

Warranty	Limited Lifetime
----------	------------------