

Bradshaw Junior Bed

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- Inspired by Drive remains a trusted, one-stop-shop for pediatric equipment and is now proud to expand its Sun Up to Sun Down collection with the addition of the Bradshaw Junior Bed. Designed specifically for children and young adults with special needs, the Bradshaw Junior combines reliable safety features with clinical-grade functionality to support everyday care routines. Delivered on a convenient transport stand and with each section weighing under 44 lb, its engineered for smooth setup while offering dependable performance with a robust 363 lb weight capacity. This fully electric bed features 360° side enclosure, handset lockbox, two-stage unlocking rails, and is tested to meet strict entrapment standards, delivering trusted protection and peace of mind for caregivers.
- Lockout handset function and dedicated safety box keep handset protected to reduce tampering
- Adjustable height and positioning support easier transfers and personalized postural alignment
- Bed frame is delivered on an integrated transport stand, making setup and relocation smooth while also doubling as a compact storage solution
- 363 lb weight capacity accommodates both children and young adults for long-term reliability as needs evolve
- All-In-One Bed Set - includes bed frame, side rails, and mattress
- Each section weighs under 44 lb and streamlined construction makes setup quick and effortless
- Reverse Trendelenburg: 0° - 12°
- Backrest Section Tilt: 0° - 74°
- Thigh Section Tilt: 0° - 31.5°
- Mattress Type: Precision-Cut Foam
- HCPCS Code: E0329
- Warranty: Limited Lifetime Frame and Welds, 5 Years Motor and Pendant, 1 Year All Other Parts

PRODUCTS	Description
BSB-ARAILPADS	Bradshaw Jr. Bed Rail Pads
BSBOAK-STD-FULL	Bradshaw Bed, Standard, Oak, w/ Lock Box
SPECIFICATIONS	Description
Warranty	Limited Lifetime Frame and Welds, 5 Years Motor and Pendant, 1 Year All Other Parts
Product Dimensions	84.6"(L) x 36.6" (W)
Weight Capacity	363 lb
Casters	6" - Dual-Locking
Color	Honey Oak
HCPCS	E0329
Motor	Full Electric w/ Hi Lo
Bed Height Minimum With Casters	36.6
Bed Height Maximum With Casters	53.1"