

FUSION™ 2 IN 1

ROLLATOR/TRANSPORT CHAIR



Two Premium Products in One for Versatile All-Day Support

The M+D® by Drive® Fusion™ delivers premium two-in-one functionality to serve more customers. The versatile design converts between rollator and transport chair in just seconds allowing for flexible use throughout the day while offering customers full-featured comfort and style.



Superior Comfort

Designed to impress discerning customers with premium features including ergonomic handgrips, a padded seat, full backrest, three footrest positions, and a removable bag that folds with the frame.



Enhanced Stability & Control

Provides enhanced stability indoors and outdoors thanks to a stable center of gravity and 8" wide, treaded wheels positioned outside the handles. Navigates bumps with ease with included curb climber.



Travels Easily

Fits seamlessly into active lifestyles with easy side-folding via pull strap and a locking frame and patented Swing n' Lock footrests that clip to the frame when folded for secure lifting and storage.



Item #	Description	UPC	Product Specs
MD700935BK	Fusion Rollator/Transport Chair, Onyx	822383030548	Handle Height 31" - 36.5"
MD700935BL	Fusion Rollator/Transport Chair, Sapphire	822383030555	Seat Height 20"
MD700935GN	Fusion Rollator/Transport Chair, Forest	822383030562	Seat Width 19.5"
			Product Width 27"
			Product Weight 21 lb

Weight Capacity	300 lb
Frame	Side-fold
Wheels	8" wide, treaded
Warranty	Limited Lifetime