




Confidence by Design

Thoughtfully engineered and beautifully designed, these bathroom safety solutions blend patented safety innovation with style that complements modern interiors.



M+D X-FOLD™ SHOWER CHAIR

Stylish, Premium Support for Your Bath



Weight
Capacity

 **Retail Packaged**

MD50100GR

Storm



MD50100WH

Pearl



MD50100TP

Stone



Stays Drier Between Uses

Built-in drainage holes prevent water pooling for a more comfortable, safer seat



Meets Your Needs

Removable backrest and arms allows for use as a chair or bench



Safer on Wet Surfaces

Anti-slip, water-dispersing feet provide strong traction to reduce sliding risk

Product Specs

Weight Capacity

300 lb

Product Weight

11.5 lb

Overall height

33.5" - 36.6"

Seat Heights

18", 20", 21"

Seat Dimensions

19.7" W x 14.5" D

Front Feet Width

17"

Rear Feet Width

15.4"

Retail Box Dimensions

17.75" x 6.6" x 24.5"

Warranty

Limited Lifetime

M+D 360™ TRANSFER BENCH

Glide In & Out with Confidence



Weight
Capacity

 **Retail Packaged**

MD50200GR Storm 

MD50200WH Pearl 



Smooth, Controlled Sliding

Large easy-glide rollers reduce friction for safer transfers



Rotates to Your Best Angle

Seat rotates with a lever and auto-locks so you can enter from your most comfortable position



Secure on Wet Surfaces

Non-slip feet grip wet floors for stability you can trust during transfers

Product Specs

Weight Capacity
Product Weight
Overall Height
Overall Length
Overall Width

350 lb
23 lb
34.25" – 40.35"
36"–38"
20"–21.7"

Seat Heights
Seat Dimensions
Bathtub Wall Height
Retail Box Dimensions
Warranty

19"- 25.5"
19.7" (W) x 13.8" (D)
14.3" – 20.6"
20" x 8.6" x 30.75"
Limited Lifetime

FOLDING SHOWER CHAIR

Stable Support That Stores Flat



 **Retail Packaged**

MD770525WH Pearl 



Perfect for Small Spaces

Folds flat for quick, compact storage in tight bathrooms, closets, or travel



Safer on Wet Surfaces

Anti-slip, water-dispersing feet provide strong traction to reduce sliding risk



Custom Height for Your Comfort

Adjustable legs from 15" to 20" help you find the most comfortable seated position

Product Specs

Weight Capacity	250 lb	Seat Dimensions	19.6" (W) x 12.2" (D)
Product Weight	6 lb	Outside Leg Dimensions	15.25" x 14.5"
Overall Height	27"-33"	Retail Box Dimensions	14" x 5.6" x 20.25"
Seat Height	15"-20"	Warranty	Limited Lifetime

SHOWER BENCH

Lightweight Strength That Adapts to You



Weight Capacity

Retail Packaged

MD770535WH Pearl



Drier, Safer Seating

Drainage holes and textured surface keep the seat dry and slip-free



Safer on Wet Surfaces

Anti-slip tips on each leg grip wet floors for dependable traction and stability



Custom Height for Your Comfort

Quick push-pin adjustment from 15" to 21"

Product Specs

Weight Capacity	300 lb	Outside Leg Dimensions	15" x 15"
Product Weight	5 lb	Retail Box Dimensions	13" x 4.25" x 22.5"
Overall Height	15"-21.5"	Warranty	Limited Lifetime
Seat Dimensions	20" (W) x 12" (D)		

MULTI TUB RAIL

Custom Support Where You Need It



**300
LBS**
Weight
Capacity

 **Retail Packaged**

MD785340WH Pearl 



Adjust to Your Perfect Grip

Rotating, lockable handgrip adjusts to 8 positions, parallel or perpendicular to the tub



Customize Your Reach

4 height adjustments ensure perfect support for every user



Fits Most Tubs

Adjusts to bathtub walls between 3" and 7" wide while protecting tub walls

Product Specs

Weight Capacity 300 lb
Product Weight 6 lb
Product Height 18" - 21"

Clamp Adjustments 3"-7"
Retail Box Dimensions 7.25" x 6.6" x 12.6"
Warranty Limited Lifetime

MASSAGE BATH MAT

Security in Every Bath



 Retail Packaged

MD785500WH Pearl 



Feel the Relief

4 invigorating massage zones designed to soothe tired feet



Stays Secure

400 suction cups grip the tub for stability you can trust



Fits Your Bath Chair

Contoured design works seamlessly with standard bath chairs

Product Specs

Weight Capacity
Product Weight
Product Dimensions

300 lb
2.8 lb
39.4" x 15.7"

Retail Box Dimensions
Warranty

4.25" x 4.25" x 16.75"
Limited Lifetime



www.drivemedical.com

©2026 Medical Depot, Inc. dba Drive Medical, 99 Seaview Boulevard, Port Washington, NY 11050. All trademarks used in association with the sale of products of Drive Medical are trademarks owned by Medical Depot, Inc. All other trademarks, trade names, service marks and logos referenced herein belong to their respective companies.

MDBATHBROCH26.VERA.4.26