

NITRO CRUISE™ Wheelchair

Contemporary Style, Premium Comfort & Intuitive Ease of Use

drive MEDICAL



K4 Design that Stands Apart

The contemporary style and premium comfort of Nitro® come to the lightweight K4 category with the Nitro® Cruise™ Wheelchair.

Built **under 30 pounds** with a **sleek frame, streamlined adjustability, terrain-ready wheels, and a unique new lumbar adjustment system**, the Cruise delivers the elevated experience active users want while remaining simple to stock, set up, and service.

**Convenient
Storage
& Transport**



Adjustability Without Complexity

Patented seat depth adjustability, height-adjustable arms, and infinitely adjustable footrests—including reversible, angle adjustable footplates—give providers an affordable, customizable chair without complex setup. The padded back cushion with removable foam lumbar and trunk support blocks adds a unique comfort upgrade customers notice immediately.

Lightweight Build That Reduces Service Burden

At just 28 pounds, the aluminum frame is easier for caregivers to lift and for users to self-propel, reducing complaints and returns related to weight.

Modern Aesthetic That Drives Demand

The sleek curved frame brings style to buyers looking beyond utilitarian styling, giving you a premium option that commands attention on the showroom floor and online.

Transport-Friendly Design

Quick-release wheels and a fold-down backrest make the Cruise easy to load, store, and recommend to customers who travel often.

Durable, Easy-to-Clean Seating System

The removable, machine-washable padded back cushion and positioning blocks simplify cleaning, making it ideal for homecare environments and long-term daily use.

Product Specifications

Skus	(see page 2)
Frame	Aluminum
Colors	Red or Black Frame
Seat Widths	16", 18", 20"
Seat Depth	16 - 18" (adjustable)
Weight Capacity	300 lb.
Product Weight	28 lb.
Axle Type	Quick-Release
Footrest Options	Swing-away Footrests or Elevating Leg Rests
Warranty	Limited lifetime on frame
HCPCS Code(s)	K0004 (pending)



Product Information

Item #	Description	Color	UPC	UOM
NW16FBADARD-SF	16", Flip-back Adjustable Desk Arm, Swing-away Footrest	●	822383032993	1/ea
NW16FBADABK-SF	16", Flip-back Adjustable Desk Arm, Swing-away Footrest	●	822383033020	1/ea
NW16FBADARD-ELR	16", Flip-back Adjustable Desk Arm, Elevating Legrests	●	822383034720	1/ea
NW16FBADABK-ELR	16", Flip-back Adjustable Desk Arm, Elevating Legrests	●	822383034751	1/ea
NW18FBADARD-SF	18", Flip-back Adjustable Desk Arm, Swing-away Footrest	●	822383033006	1/ea
NW18FBADABK-SF	18", Flip-back Adjustable Desk Arm, Swing-away Footrest	●	822383033037	1/ea
NW18FBADARD-ELR	18", Flip-back Adjustable Desk Arm, Elevating Legrests	●	822383034737	1/ea
NW18FBADABK-ELR	18", Flip-back Adjustable Desk Arm, Elevating Legrests	●	822383034768	1/ea
NW20FBADARD-SF	20", Flip-back Adjustable Desk Arm, Swing-away Footrest	●	822383033013	1/ea
NW20FBADABK-SF	20", Flip-back Adjustable Desk Arm, Swing-away Footrest	●	822383033044	1/ea
NW20FBADARD-ELR	20", Flip-back Adjustable Desk Arm, Elevating Legrests	●	822383034744	1/ea
NW20FBADABK-ELR	20", Flip-back Adjustable Desk Arm, Elevating Legrests	●	822383034775	1/ea

Accessories



Auto-Clasp Seatbelt
Item# STDS850



Wheel Lock Extension
Item# STDS801



Cup holder
Item# 102662-CUP



DualEase Gel Cushion
Item# 14991, 14992, 14993