

Nitro HD Rollator

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D




FIGURE E



FIGURE F

- Attractive, Euro-style design
- Brake cable inside frame for added safety
- Handle height easily adjusts with unique push button
- Unique suspension system for support and comfort
- Back support height easily adjusts with tool-free thumb screw
- Caster fork design enhances turning radius
- Large 10" front casters allow optimal steering and rolling comfort
- Easily folds with one hand to ultra-compact size for storage
- Under seat locking mechanism holds the Nitro closed when frame is folded
- Cross-brace design allows for side-to-side folding and added stability
- Seat is durable and comfortable
- Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded
- Lightweight, aluminum frame

PRODUCTS	Description
10266HD-BK	HD Nitro Rollator, Black, 1/ea
SPECIFICATIONS	Description
Warranty	Frame: Limited Lifetime Brake Handle and Housing: 5 Years
Seat Dimensions	23" (W) x 10" (D)
Handle Height	36" - 41"
Weight Capacity	450 lb
Seat To Floor Height	23.5"
Weight	20.5 lb
ACCESSORIES	Description
10266-BRAKE ENH	Nitro Slow Down Brake, ENH, 1PR, 25/CS, EA