

Phoenix LT 4-Wheel Scooter

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

The Phoenix 4-wheel travel scooter combines a 350 lbs weight capacity with the take-apart portability of a travel scooter. The 4-wheel design and anti-tip wheels assist with overall stability for a more secure experience. This scooter also has a height-adjustable seat with flip-back and adjustable-width arms to accommodate varying body types. The seat also helps with portability by folding flat and allowing for easier transporting or storage. Your Phoenix 4-Wheel scooter can easily and quickly disassemble into five pieces.

- Interchangeable color panels
- Seat swivels 180 degrees for easy transfer
- Headlight for increased visibility
- Breaks down to 5 pieces for easy transport
- Comes with rugged 9"x3" wheels
- Padded Delta tiller

PRODUCTS	Description
PHOENIXLT4	Phoenix LT 4-Wheel Travel POV 18" Seat
PHOENIXLT4-20	Phoenix LT 4-Wheel Travel POV 20" Seat
SPECIFICATIONS	Description
Warranty	24 Month Limited Warranty on Electronics
Product Weight	127 lb
Product Dimensions	42.3" (L) x 20.1" (W) x 36.6" (H)
Weight Capacity	350 lb
Top Speed	4 mph
Batteries (pair)	12V x 22AH
Battery Charger	2A Offboard
Climbing Angle	9°
Color	Gray, Red, Blue
Controller	PG 45A
Ground Clearance	3.9"
HCPCS	K0801
Maximum Range	14 miles
Motor	24V x 250W x 170
Seat Type	Stadium
Turning Radius	70"
Warranty on Batteries	12 Month Warranty on Batteries
Warranty on Frame	Limited Lifetime Warranty on Frame

SPECIFICATIONS**Description**

Wheel Size	Front: 9" x 3" Rear: 9" x 3"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]