

## Panther 4-Wheel Heavy Duty Scooter

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



The Panther scooter offers modern styling and a host of excellent features resulting in a smooth and comfortable driving experience. Just some of these great features include one touch infinite tiller adjustment and delta controls for easy steering and throttle control, a plush Captains Seat that is height adjustable and features a reclining backrest with adjustable headrest, 14" alloy wheels on front and rear suspension for a smooth ride on all terrain, extra bright LED headlights and turn signals, dual rear view mirrors, and a large front basket with carry handle. For the ultimate driving experience dont settle for anything less than the Panther.

- Front and rear suspension on 14" pneumatic tires for the smoothest ride possible
- One-touch infinite tiller adjustment
- Built-in USB port for charging your mobile devices
- Extra-bright LED lighting for superb visibility
- Ergonomic delta control handles for ease of steering and throttle control
- Large front plastic basket with carry handle
- Captain's seat is height adjustable and features a fold-down, reclining backrest with adjustable headrest
- LED headlights, parking lights, and turn signals for superior visibility
- Dual rear view mirrors
- Anti-tip wheels
- Easy freewheel mode

| PRODUCTS           | Description                              |
|--------------------|------------------------------------------|
| PANTHER22CS        | Panther HD Scooter, 22" Capt.Seat,1/ea   |
| SPECIFICATIONS     | Description                              |
| Warranty           | 24 Month Limited Warranty on Electronics |
| Product Weight     | 264 lb                                   |
| Product Dimensions | 56" (L) x 26" (W) x 53" (H)              |
| Seat Dimensions    | 20" (D) x 22" (W)                        |
| Weight Capacity    | 425 lb                                   |
| Top Speed          | 8 mph                                    |
| Batteries (pair)   | 12V x 50AH                               |
| Battery Charger    | 5A Offboard                              |
| Climbing Angle     | 7.5°                                     |
| Color              | Red                                      |
| Controller         | Dynamic Rhino 2 120A                     |
| Ground Clearance   | 4"                                       |
| HCPCS              | N/A                                      |
| Maximum Range      | 25 miles                                 |

| SPECIFICATIONS        | Description                        |
|-----------------------|------------------------------------|
| Motor                 | 24V x 800W x 4500 RPM              |
| Seat Type             | Captain - Reclining High           |
| Turning Radius        | 53"                                |
| Warranty on Batteries | 12 Month Warranty on Batteries     |
| Warranty on Frame     | Limited Lifetime Warranty on Frame |
| Wheel Size            | Front: 14" x 3"   Rear: 14" x 3"   |