

Prowler 4-Wheel Mobility Scooter

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Delta tiller with ergonomic throttle control for better grip and additional steering comfort
- Height-adjustable swivel seat features fold-down backrest making storage easier.
- Seat back reclines and features an adjustable headrest
- Seat position can be adjusted forward or backward to maximize comfort specific to each user.
- Pneumatic tires and full suspension allow for a smooth ride
- Large plastic carry basket to take items on-the-go
- Easy-to-adjust tiller
- Programmable controller to adjust driving parameters based on user's needs
- Anti-tip wheels for additional safety
- Front headlight offers increased visibility
- Rearview mirrors
- Easy Freewheel operation
- Additional hand-operated drum brake

PRODUCTS

Description

PROWLER3410MG22CS

PROWLER 3410 BLACK 22" CAPT SEAT