

Aluminum Folding Canes with Gel Grip, Height Adjustable

HCPCS: E0100



FIGURE A




FIGURE B



FIGURE C



FIGURE D



FIGURE E

- New gel grip reduces hand stress and fatigue
- Aluminum folding cane, height adjustable
- Cane folds into four convenient parts for easy storage
- Manufactured with sturdy, extruded aluminum tubing
- Comes with plastic clip to hang cane when folded
- Available in Red Crackle, Blue Crackle, Black, and Hearts
- Handle height adjusts from 32" to 36"
- Weight Capacity 300 lb

PRODUCTS

Description

RTL10370BC	Cane,Folding,Gel Grip,BLCrackle,1EA,6/CS
RTL10370BK	Cane,Folding,Gel Grip,BLK,1EA,6/CS
RTL10370RC	Cane,Folding,Gel Grip,Red Crkle,1EA,6/CS

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Handle Height	32" - 36"
Weight Capacity	300 lb